

3-PLY SURGICAL MASK



DESCRIPTION

- Pleated 3 layer design, the filter efficiency is higher than 99%
- Soft lined, water resistant, non-woven material
- Latex and fiberglass free
- Adjustable nose piece and ear-loop design for convenience and comfort

3-PLY SURGICAL MASK



SPECIFICATIONS

Height	6.9" (17.5cm)
Width	3.7" (9.5cm)
Colour	Blue

OTHER INFORMATION

Material	Polypropylene Non-woven fabric
Filter Rate	BFE ≥ 99%
Use Range	Hospitals/medical (Personal protection)
Disposable	Yes
Type	Sterility (Surgical) Non-sterility (Medical)
Style	Ear-loop

