

## 3-PLY SURGICAL MASK (KIDS)



### DESCRIPTION

- Pleated 3 layer design, the filter efficiency is higher than 99%
- Soft lined, water resistant, non-woven material
- Latex and fiberglass free
- Adjustable nose piece and ear-loop design for convenience and comfort
- Fits most children (Available in blue or pink)

# 3-PLY SURGICAL MASK (KIDS)

## SPECIFICATIONS



Height	6.5" (15.24cm)
Width	3.5" (8.89cm)
Colour	Blue or Pink

## OTHER INFORMATION

Material	Polypropylene Non-woven fabric
Filter Rate	BFE ≥ 99%
Use Range	Hospitals/medical (Personal protection)
Disposable	Yes
Type	Sterility (Surgical) Non-sterility (Medical)
Style	Ear-loop



# 3-PLY SURGICAL MASK (KIDS)

WORLD VISION INITIATIVE

*With the purchase of a pack of 5 kids masks, Dent-X Canada will donate 3 masks to World Vision, to help keep those less fortunate both safe and healthy.*



**DENT-X**  
CANADA  
1-866-791-0909