


## Procedure Masks 3-Ply Kids Level 3



### Description & Classification

Model	3-Ply Kids, Level 3
Product Name	Procedure Mask
Strap Type	Earloop
Size	Standard, Fits most kids
Sample	
Measurements	W: 5.7" X H: 3.5"
Colour	Blue/Pink
Manufactured In	Canada
Use Range	Medical Facilities Personal Protection

### Material and Performance Specification

Material	Non-woven Fabric
Contains Latex	No
Bacterial Filter Efficiency	≥99.78%
Differential Pressure	≤3.78 mmH <sub>2</sub> O
Particulate Filter Efficiency	≥99.96%
Fluid Resistance	≥160mm Hg
Flammability Rating	Class I

### Performance Requirements

Test	ASTM F2100			3-Ply Kids, Level 3
	Level 1	Level 2	Level 3	
Bacterial Filtration Efficiency	≥95%	≥98%	≥98%	Pass
Differential Pressure	< 5 mmH <sub>2</sub> O	< 6 mmH <sub>2</sub> O	< 6 mmH <sub>2</sub> O	Pass
Particulate Filtration Efficiency	≥95%	≥98%	≥98%	Pass
Fluid Resistance	80 mm Hg	120 mm Hg	160 mm Hg	Pass
Flammability Rating	Class 1	Class 1	Class 1	Pass, Class 1



### Packaging Specification

Mask Weight	0.11 oz	0.11 oz
Masks in Box	10	50
Weight of Box with Masks	0.11 lbs	0.38 lbs
Box Dimensions	L: 8.25" x W: 5" x H: 1"	L: 7.4" x W: 3.9" x H: 3.9"
Masks Per Case	500	1000
Boxes Per Case	50	20
Weight of Case	6.5 lbs	8.6 lbs
Case Dimensions	L: 20" x W: 15.5" x H: 8.75"	L: 20" x W: 15.5" x H: 8.75"
Masks Per Pallet	21000	42000
Cases Per Pallet	42	42
Pallet Dimensions	L: 41" x W: 48" x H: 66"	L: 41" x W: 48" x H: 66"